

# ERG 910 PLUS ERG 911 PLUS

Functional, high-quality ergometers



Tolmed



**SCHILLER**  
The Art of Diagnostics

# ERG 910 PLUS

# ERG 911 PLUS

- ⌘ Pivoting display 180°
- ⌘ Easy to operate
- ⌘ Almost noiseless operation
- ⌘ Low access
- ⌘ Stable steel tube construction



## APPLICATION AREAS

The bike ergometers ERG 910 plus and ERG 911 plus were developed for use in the cardiology sector based on the proven previous models, with particular care to ensure highest stability and practicability. The devices are used for performing defined exercises during ergometric examinations and therapeutic applications.

ERG 911 plus is available in both versions "easy" and "smart". "Easy" devices are intended to be used with a master device (remote operation: ECG, PC, etc.) whereas "smart" models can be used as stand-alone devices. Furthermore, many options can be added to the available features on both devices.

Regardless of their features, the devices meet the highest quality standards in compliance with DIN VDE 0750-238 for accurate physical exertion tests to conduct measurements in cardiovascular and pulmonary function diagnostics.

- ⌘ Blood pressure measurement using QRS trigger
- ⌘ Electrical saddle height adjustment
- ⌘ SpO<sub>2</sub> measurement
- ⌘ Bluetooth connection
- ⌘ ECG amplifier module



### ERG 910 PLUS - EASY

- ❖ Remote operation (type "easy")
- ❖ load range 20-500 Watts
- ❖ optional: power supply 110-120V
- ❖ optional: 20-999 Watt



### ERG 911 PLUS - EASY / SMART

- ❖ remote (type "easy") and stand-alone (type "smart") operation
- ❖ load range 20-999 Watts
- ❖ various optional extensions available for both models (see technical data)



### EASY

Ergometers equipped with an "easy" control panel, on which they receive exercise specifications from a master device (ECG, PC, ...), are used for remote operation. All major exercise data (such as blood pressure, heart rate, pedalling rate, instruction concerning correct pedalling rate) is displayed in a compact format. The optional features, such as blood pressure measurement and electric saddle height adjustment, can also be conveniently adjusted from the Touch display.



### SMART

Ergometers equipped with a "smart" control panel, can be used for remote as well as for stand-alone operation. Exercise instructions can be transmitted from a master device (ECG, PC, ...) or from a programmable ergometry program located on the ergometer. Further, training programs with exercise types such as pulse steady state, constant load and interval training can also be programmed.



**Americas**  
SCHILLER Americas Inc.  
Doral, Florida 33172

**North America:**

Phone +1 786 845 06 20  
Fax +1 786 845 06 02  
sales@schilleramerica.com  
www.schilleramerica.com

**Latin America & Caribbean**

Phone +1 305 591 11 21  
Fax +1 786 845 06 02  
sales@schillerla.ch  
www.schillerla.ch



**Asia**  
SCHILLER Asia-Pacific / Malaysia  
52200 Kuala Lumpur, Malaysia  
Phone +603 6272 3033  
Fax +603 6272 2030  
sales@schiller.com.my  
www.schiller-asia.com



**Austria**  
SCHILLER Handelsgesellschaft m.b.H.  
A-4040 Linz  
Phone +43 732 70 99  
Fax +43 732 757 000  
sales@schiller.at  
www.schiller.at



**China**  
Alfred Schiller (Beijing) Medical Equipment  
Co.,Ltd.  
100015 Beijing, China  
Phone +86 010 52007020  
Fax +86 010 52007020-8016  
info@schillerchina.com  
www.schillermedical.cn



**Croatia**  
SCHILLER d.o.o.  
10000 Zagreb  
Phone +385 1 309 66 59  
Fax +385 1 309 66 60  
info@schillerzg.hr  
www.schiller.hr



**France**  
SCHILLER Médical  
F-67160 Wissembourg  
Phone +33 3 88 63 36 00  
Fax +33 3 88 63 36 49  
info@schiller.fr  
www.schiller-medical.com



**France (distribution France)**  
SCHILLER France S.A.S.  
F-77608 Bussy St Georges  
Phone +33 1 64 66 50 00  
Fax +33 1 64 66 50 10  
contact@schillerfrance.fr  
www.schiller-france.com



**Germany**  
SCHILLER Medizintechnik GmbH  
D-85622 Feldkirchen b. München  
Phone +49 89 62 99 81 0  
Fax +49 89 609 50 90  
info@schillermed.de  
www.schillermed.de



**Hungary**  
SCHILLER Diamed Ltd.  
H-1141 Budapest  
Phone +36 1 383 4780 / +36 1 460 9491  
Fax +36 1 383 4778  
info@schillerhungary.hu  
www.schillerhungary.hu



**India**  
SCHILLER Healthcare India Pvt. Ltd.  
Mumbai – 400 059, India  
Phone +91 22 6152 3333 / 2920 9141  
Fax +91 22 2920 9142  
sales@schillerindia.com  
www.schillerindia.com



**Poland**  
SCHILLER Poland Sp. z o.o.  
PL-02-729 Warszawa  
Phone +48 22 843 20 89 / +48 22 647 35 90  
Fax +48 22 843 20 89  
schiller@schiller.pl  
www.schiller.pl



**Russia & C.I.S.**  
AO SCHILLER.RU  
125040 Moscow, Russia  
Phone +7 (495) 970 11 33  
Fax +7 (495) 956 29 10  
mail@schiller.ru  
www.schiller.ru



**Serbia**  
SCHILLER d.o.o.  
11010 Beograd  
Phone +381 11 39 79 508  
Fax +381 11 39 79 518  
info@schiller.rs  
www.schiller.rs



**Slovenia**  
SCHILLER d.o.o.  
2310 Slovenska Bistrica  
Phone +386 2 843 00 56  
Fax +386 2 843 00 57  
info@schiller.si  
www.schiller.si



**Spain**  
SCHILLER ESPAÑA, S.A.  
ES-28232 Las Rozas/Madrid  
Phone +34 91 713 01 76  
Fax +34 91 355 79 33  
schiller@schiller.es  
www.schiller.es



**Switzerland**  
SCHILLER-Reomed AG  
CH-8953 Dietikon  
Phone +41 44 744 30 00  
Fax +41 44 740 37 10  
sales@schiller-reomed.ch  
www.schiller-reomed.ch



**Turkey**  
SCHILLER TÜRKİYE  
Okmeydanı-Sisli – Istanbul  
Phone +90 212 210 8681 (pbx)  
Fax +90 212 210 8684  
info@schiller.com.tr  
www.schiller-turkiye.com



**United Kingdom**  
SCHILLER UK  
Bellshill, ML4 3PR  
Phone +44 1698 744 505  
Fax +44 1698 744 474  
sales@schilleruk.com  
www.schilleruk.com

**Van der Tol medische apparaten v.o.f.**  
Tasveld 54 8271 RW IJsselmuiden  
Tel.: 038 333 83 80 Fax: 038 333 86 40  
[www.tolmed.nl](http://www.tolmed.nl) [info@tolmed.nl](mailto:info@tolmed.nl)

Article no.: 2.500772  
2018/03\_EN



**Tolmed**



**SCHILLER**

The Art of Diagnostics