

CARDIOLINE

Every single beat.



Walk Free

Wearable 3 channel cable-free ECG
Holter patch recorder

- Cable-free 3 channel holter (Leads I, II, III)
- Discrete, Small & Compact as well as lightweight: less than 30g
- Standard Electrode Application
- Bluetooth Realtime ECG Waveform Hook Up / Quality Check
- Recording up to 9 days with a single battery





Walk Free

Technical specifications

Leads	3 Leads: I, II, III
ADC Resolution	16 bit / 2,5 mV LSB
Sampling Rate	250 samples/sec
Frequency Bandwidth	0.05 – 70 Hz
Recording Duration	Up to 9 days with a single battery
Data transfer	Bluetooth (4.0) and USB
Regulations	93/42/EEC, EN 60601-1: Medical electrical equipment; defibrillation protected
CE Product Classification	Class IIa
IP rating	IP54B
Power	AAAA single battery
Dimensions	71,5 x 46 x 14,7 mm
Weight	30 g

